

Antibodies Product Datasheet

Product Name:	Galectin-4 Antibody
Catalog Number:	C13077
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LEG4. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	LGALS4
Protein Name:	Galectin-4
Other Name:	LGALS4; Galectin-4; Gal-4; Antigen NY-CO-27; L-36 lactose-binding protein; L36LBP; Lactose-binding lectin 4
Human Swiss Prot No.:	P56470
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Galectin-4 Polyclonal Antibody detects endogenous levels of Galectin-4 protein.
Storage and Stability:	-20°C/1 year