

Antibodies Product Datasheet

Product Name:	Galectin-8 Antibody
Catalog Number:	C13079
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LEG8. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	LGALS8
Protein Name:	Galectin-8
Other Name:	LGALS8; Galectin-8; Gal-8; Po66 carbohydrate-binding protein; Po66-CBP; Prostate carcinoma tumor antigen 1; PCTA-1
Human Swiss Prot No.:	O00214
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Galectin-8 Polyclonal Antibody detects endogenous levels of Galectin-8 protein.
Storage and Stability:	-20°C/1 year