

Antibodies Product Datasheet

Product Name:	MIP-T3 Antibody
Catalog Number:	C13085
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MIPT3. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	TRAF3IP1
Protein Name:	TRAF3-interacting protein 1
Other Name:	TRAF3IP1; MIPT3; TRAF3-interacting protein 1; Interleukin-13 receptor alpha 1-binding protein 1; Microtubule-interacting protein associated with TRAF3; MIP-T3
Human Swiss Prot No.:	Q8TDR0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MIP-T3 Polyclonal Antibody detects endogenous levels of MIP-T3 protein.
Storage and Stability:	-20°C/1 year