

Antibodies Product Datasheet

Product Name:	MSH3 Antibody
Catalog Number:	C13090
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MSH3. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	MSH3
Protein Name:	DNA mismatch repair protein Msh3
Other Name:	MSH3; DUC1; DUG; DNA mismatch repair protein Msh3; hMSH3; Divergent upstream protein; DUP; Mismatch repair protein 1; MRP1
Human Swiss Prot No.:	P20585
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MSH3 Polyclonal Antibody detects endogenous levels of MSH3 protein.
Storage and Stability:	-20°C/1 year