

Antibodies Product Datasheet

Product Name:	Rad23A Antibody
Catalog Number:	C13113
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RAD23A. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	RAD23A
Protein Name:	UV excision repair protein RAD23 homolog A
Other Name:	RAD23A; UV excision repair protein RAD23 homolog A; HR23A; hHR23A
Human Swiss Prot No.:	P54725
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Rad23A Polyclonal Antibody detects endogenous levels of Rad23A protein.
Storage and Stability:	-20°C/1 year