

Antibodies Product Datasheet

Product Name:	KLF7 Antibody
Catalog Number:	C13415
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KLF7. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	KLF7
Protein Name:	Krueppel-like factor 7
Other Name:	KLF7; UKLF; Krueppel-like factor 7; Ubiquitous krueppel-like factor
Human Swiss Prot No.:	O75840
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KLF7 Polyclonal Antibody detects endogenous levels of KLF7 protein.
Storage and Stability:	-20°C/1 year