

Antibodies Product Datasheet

Product Name:	Maf1 Antibody
Catalog Number:	C13419
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Maf1 . at AA range: 110-190
Concentration:	1 mg/ml
Gene Name:	MAF1
Protein Name:	Repressor of RNA polymerase III transcription MAF1 homolog
Other Name:	MAF1; Repressor of RNA polymerase III transcription MAF1 homolog
Human Swiss Prot No.:	Q9H063
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Maf1 Polyclonal Antibody detects endogenous levels of Maf1 protein.
Storage and Stability:	-20°C/1 year