

Antibodies Product Datasheet

Product Name:	LPAAT- δ Antibody
Catalog Number:	C14005
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AGPAT4. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	AGPAT4
Protein Name:	1-acyl-sn-glycerol-3-phosphate acyltransferase delta
Other Name:	AGPAT4; 1-acyl-sn-glycerol-3-phosphate acyltransferase delta; 1-acylglycerol-3-phosphate O-acyltransferase 4; 1-AGP acyltransferase 4; 1-AGPAT 4; Lysophosphatidic acid acyltransferase delta; LPAAT-delta
Human Swiss Prot No.:	Q9NRZ5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LPAAT- δ Polyclonal Antibody detects endogenous levels of LPAAT- δ protein.
Storage and Stability:	-20°C/1 year