

Antibodies Product Datasheet

Product Name:	LPAAT- ϵ Antibody
Catalog Number:	C14006
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AGPAT5. AA range:241-290
Concentration:	1 mg/ml
Gene Name:	AGPAT5
Protein Name:	1-acyl-sn-glycerol-3-phosphate acyltransferase epsilon
Other Name:	AGPAT5; 1-acyl-sn-glycerol-3-phosphate acyltransferase epsilon; 1-acylglycerol-3-phosphate O-acyltransferase 5; 1-AGP acyltransferase 5; 1-AGPAT 5; Lysophosphatidic acid acyltransferase epsilon; LPAAT-epsilon
Human Swiss Prot No.:	Q9NUQ2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LPAAT- ϵ Polyclonal Antibody detects endogenous levels of LPAAT- ϵ protein.
Storage and Stability:	-20°C/1 year