

Antibodies Product Datasheet

Product Name:	LPAAT- γ Antibody
Catalog Number:	C14008
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AGPAT3. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	AGPAT3
Protein Name:	1-acyl-sn-glycerol-3-phosphate acyltransferase gamma
Other Name:	AGPAT3; 1-acyl-sn-glycerol-3-phosphate acyltransferase gamma; 1-acylglycerol-3-phosphate O-acyltransferase 3; 1-AGP acyltransferase 3; 1-AGPAT 3; Lysophosphatidic acid acyltransferase gamma; LPAAT-gamma
Human Swiss Prot No.:	Q9NRZ7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LPAAT- γ Polyclonal Antibody detects endogenous levels of LPAAT- γ protein.
Storage and Stability:	-20°C/1 year