

Antibodies Product Datasheet

Product Name:	MRP-S7 Antibody
Catalog Number:	C14044
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRPS7. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	MRPS7
Protein Name:	28S ribosomal protein S7 mitochondrial
Other Name:	MRPS7; 28S ribosomal protein S7; mitochondrial; MRP-S7; S7mt; bMRP-27a; bMRP27a
Human Swiss Prot No.:	Q9Y2R9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRP-S7 Polyclonal Antibody detects endogenous levels of MRP-S7 protein.
Storage and Stability:	-20°C/1 year