

Antibodies Product Datasheet

Product Name:	3 β -HSD7 Antibody
Catalog Number:	C14049
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HSD3B7. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	HSD3B7
Protein Name:	3 beta-hydroxysteroid dehydrogenase type 7
Other Name:	HSD3B7; 3 beta-hydroxysteroid dehydrogenase type 7; 3 beta-hydroxysteroid dehydrogenase type VII; 3-beta-HSD VII; 3-beta-hydroxy-Delta(5)-C27 steroid oxidoreductase; C(27) 3-beta-HSD; Cholest-5-ene-3-beta; 7-alpha-diol 3-beta-dehydrogenase
Human Swiss Prot No.:	Q9H2F3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	3 β -HSD7 Polyclonal Antibody detects endogenous levels of 3 β -HSD7 protein.
Storage and Stability:	-20°C/1 year