

Antibodies Product Datasheet

Product Name:	MRP-L34 Antibody
Catalog Number:	C14072
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRPL34. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	MRPL34
Protein Name:	39S ribosomal protein L34 mitochondrial
Other Name:	MRPL34; 39S ribosomal protein L34; mitochondrial; L34mt; MRP-L34
Human Swiss Prot No.:	Q9BQ48
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRP-L34 Polyclonal Antibody detects endogenous levels of MRP-L34 protein.
Storage and Stability:	-20°C/1 year