

Antibodies Product Datasheet

Product Name:	MRP-L40 Antibody
Catalog Number:	C14077
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRPL40. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	MRPL40
Protein Name:	39S ribosomal protein L40 mitochondrial
Other Name:	MRPL40; NLVCF; URIM; 39S ribosomal protein L40; mitochondrial; L40mt; MRP-L40; Nuclear localization signal-containing protein deleted in velocardiofacial syndrome; Up-regulated in metastasis
Human Swiss Prot No.:	Q9NQ50
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRP-L40 Polyclonal Antibody detects endogenous levels of MRP-L40 protein.
Storage and Stability:	-20°C/1 year