

Antibodies Product Datasheet

Product Name:	MRP-L47 Antibody
Catalog Number:	C14083
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRPL47. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	MRPL47
Protein Name:	39S ribosomal protein L47 mitochondrial
Other Name:	MRPL47; NCM1; CGI-204; 39S ribosomal protein L47; mitochondrial; L47mt; MRP-L47; Nasopharyngeal carcinoma metastasis-related protein 1
Human Swiss Prot No.:	Q9HD33
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRP-L47 Polyclonal Antibody detects endogenous levels of MRP-L47 protein.
Storage and Stability:	-20°C/1 year