

Antibodies Product Datasheet

Product Name:	HIBADH Antibody
Catalog Number:	C14092
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HIBADH. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	HIBADH
Protein Name:	3-hydroxyisobutyrate dehydrogenase mitochondrial
Other Name:	HIBADH; 3-hydroxyisobutyrate dehydrogenase; mitochondrial; HIBADH
Human Swiss Prot No.:	P31937
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HIBADH Polyclonal Antibody detects endogenous levels of HIBADH protein.
Storage and Stability:	-20°C/1 year