

Antibodies Product Datasheet

Product Name:	AZI1 Antibody
Catalog Number:	C14137
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AZI1. AA range:311-360
Concentration:	1 mg/ml
Gene Name:	AZI1
Protein Name:	5-azacytidine-induced protein 1
Other Name:	AZI1; CEP131; KIAA1118; 5-azacytidine-induced protein 1; Centrosomal protein of 131 kDa; Cep131; Pre-acrosome localization protein 1
Human Swiss Prot No.:	Q9UPN4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AZI1 Polyclonal Antibody detects endogenous levels of AZI1 protein.
Storage and Stability:	-20°C/1 year