

Antibodies Product Datasheet

Product Name:	RO60 Antibody
Catalog Number:	C14140
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TROVE2. AA range:311-360
Concentration:	1 mg/ml
Gene Name:	TROVE2
Protein Name:	60 kDa SS-A/Ro ribonucleoprotein
Other Name:	TROVE2; RO60; SSA2; 60 kDa SS-A/Ro ribonucleoprotein; 60 kDa Ro protein; 60 kDa ribonucleoprotein Ro; RoRNP; Ro 60 kDa autoantigen; Sjogren syndrome antigen A2; Sjogren syndrome type A antigen; SS-A; TROVE domain family member 2
Human Swiss Prot No.:	P10155
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RO60 Polyclonal Antibody detects endogenous levels of RO60 protein.
Storage and Stability:	-20°C/1 year