

Antibodies Product Datasheet

Product Name:	Ribosomal Protein L14 Antibody
Catalog Number:	C14151
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RPL14. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	RPL14
Protein Name:	60S ribosomal protein L14
Other Name:	RPL14; 60S ribosomal protein L14; CAG-ISL 7
Human Swiss Prot No.:	P50914
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ribosomal Protein L14 Polyclonal Antibody detects endogenous levels of Ribosomal Protein L14 protein.
Storage and Stability:	-20°C/1 year