

Antibodies Product Datasheet

Product Name:	Ribosomal Protein L37A Antibody
Catalog Number:	C14179
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RPL37A. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	RPL37A
Protein Name:	60S ribosomal protein L37a
Other Name:	RPL37A; 60S ribosomal protein L37a
Human Swiss Prot No.:	P61513
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ribosomal Protein L37A Polyclonal Antibody detects endogenous levels of Ribosomal Protein L37A protein.
Storage and Stability:	-20°C/1 year