

Antibodies Product Datasheet

Product Name:	ABHD11 Antibody
Catalog Number:	C14208
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ABHD11. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	ABHD11
Protein Name:	Abhydrolase domain-containing protein 11
Other Name:	ABHD11; WBSCR21; PP1226; Abhydrolase domain-containing protein 11; Williams-Beuren syndrome chromosomal region 21 protein
Human Swiss Prot No.:	Q8NFV4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ABHD11 Polyclonal Antibody detects endogenous levels of ABHD11 protein.
Storage and Stability:	-20°C/1 year