

Antibodies Product Datasheet

Product Name:	ABHD14A Antibody
Catalog Number:	C14212
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ABHD14A. AA range:201-250
Concentration:	1 mg/ml
Gene Name:	ABHD14A
Protein Name:	Abhydrolase domain-containing protein 14A
Other Name:	ABHD14A; Abhydrolase domain-containing protein 14A
Human Swiss Prot No.:	Q9BUJ0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ABHD14A Polyclonal Antibody detects endogenous levels of ABHD14A protein.
Storage and Stability:	-20°C/1 year