

Antibodies Product Datasheet

Product Name:	ABHD14B Antibody
Catalog Number:	C14213
Reactivity:	Human, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ABHD14B. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	ABHD14B
Protein Name:	Abhydrolase domain-containing protein 14B
Other Name:	ABHD14B; CIB; Abhydrolase domain-containing protein 14B; CCG1-interacting factor B
Human Swiss Prot No.:	Q96IU4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ABHD14B Polyclonal Antibody detects endogenous levels of ABHD14B protein.
Storage and Stability:	-20°C/1 year