

Antibodies Product Datasheet

Product Name:	ABHD2 Antibody
Catalog Number:	C14214
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ABHD2. AA range:251-300
Concentration:	1 mg/ml
Gene Name:	ABHD2
Protein Name:	Abhydrolase domain-containing protein 2
Other Name:	ABHD2; LABH2; Abhydrolase domain-containing protein 2; Lung alpha/beta hydrolase 2; Protein PHPS1-2
Human Swiss Prot No.:	P08910
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ABHD2 Polyclonal Antibody detects endogenous levels of ABHD2 protein.
Storage and Stability:	-20°C/1 year