

Antibodies Product Datasheet

Product Name:	PHAPI2 Antibody
Catalog Number:	C14238
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ANP32B. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	ANP32B
Protein Name:	Acidic leucine-rich nuclear phosphoprotein 32 family member B
Other Name:	ANP32B; APRIL; PHAPI2; Acidic leucine-rich nuclear phosphoprotein 32 family member B; Acidic protein rich in leucines; Putative HLA-DR-associated protein I-2; PHAPI2; Silver-stainable protein SSP29
Human Swiss Prot No.:	Q92688
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PHAPI2 Polyclonal Antibody detects endogenous levels of PHAPI2 protein.
Storage and Stability:	-20°C/1 year