

Antibodies Product Datasheet

Product Name:	ACOT12 Antibody
Catalog Number:	C14270
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ACOT12. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	ACOT12
Protein Name:	Acyl-coenzyme A thioesterase 12
Other Name:	ACOT12; CACH; CACH1; STARD15; Acyl-coenzyme A thioesterase 12; Acyl-CoA thioesterase 12; Acyl-CoA thioester hydrolase 12; Cytoplasmic acetyl-CoA hydrolase 1; CACH-1; hCACH-1; START domain-containing protein 15; StARD15
Human Swiss Prot No.:	Q8WYK0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ACOT12 Polyclonal Antibody detects endogenous levels of ACOT12 protein.
Storage and Stability:	-20°C/1 year