

Antibodies Product Datasheet

Product Name:	ACOT2 Antibody
Catalog Number:	C14271
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ACOT2. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	ACOT2
Protein Name:	Acyl-coenzyme A thioesterase 2 mitochondrial
Other Name:	ACOT2; PTE2; PTE2A; Acyl-coenzyme A thioesterase 2; mitochondrial; Acyl-CoA thioesterase 2; Acyl-coenzyme A thioester hydrolase 2a; CTE-1a; Long-chain acyl-CoA thioesterase 2; ZAP128
Human Swiss Prot No.:	P49753
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ACOT2 Polyclonal Antibody detects endogenous levels of ACOT2 protein.
Storage and Stability:	-20°C/1 year