

Antibodies Product Datasheet

Product Name:	ACOT9 Antibody
Catalog Number:	C14274
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ACOT9 . at AA range: 240-320
Concentration:	1 mg/ml
Gene Name:	ACOT9
Protein Name:	Acyl-coenzyme A thioesterase 9 mitochondrial
Other Name:	ACOT9; CGI-16; Acyl-coenzyme A thioesterase 9; mitochondrial; Acyl-CoA thioesterase 9; Acyl-CoA thioester hydrolase 9
Human Swiss Prot No.:	Q9Y305
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ACOT9 Polyclonal Antibody detects endogenous levels of ACOT9 protein.
Storage and Stability:	-20°C/1 year