

Antibodies Product Datasheet

Product Name:	ACYP2 Antibody
Catalog Number:	C14278
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ACYP2 . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	ACYP2
Protein Name:	Acylphosphatase-2
Other Name:	ACYP2; ACYP; Acylphosphatase-2; Acylphosphatase; muscle type isozyme; Acylphosphate phosphohydrolase 2
Human Swiss Prot No.:	P14621
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ACYP2 Polyclonal Antibody detects endogenous levels of ACYP2 protein.
Storage and Stability:	-20°C/1 year