

Antibodies Product Datasheet

Product Name:	LPCAT2 Antibody
Catalog Number:	C14281
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LPCAT2. AA range:321-370
Concentration:	1 mg/ml
Gene Name:	LPCAT2
Protein Name:	Lysophosphatidylcholine acyltransferase 2
Other Name:	LPCAT2; AYTL1; Lysophosphatidylcholine acyltransferase 2; LPC acyltransferase 2; LPCAT-2; LysoPC acyltransferase 2; 1-acylglycerophosphocholine O-acyltransferase; 1-alkylglycerophosphocholine O-acetyltransferase; Acetyl-CoA:lyso-platelet-ac
Human Swiss Prot No.:	Q7L5N7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LPCAT2 Polyclonal Antibody detects endogenous levels of LPCAT2 protein.
Storage and Stability:	-20°C/1 year