

Antibodies Product Datasheet

Product Name:	ADAMTS-1 Antibody
Catalog Number:	C14302
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ADAMTS-1 . at AA range: 160-240
Concentration:	1 mg/ml
Gene Name:	ADAMTS1
Protein Name:	A disintegrin and metalloproteinase with thrombospondin motifs 1
Other Name:	ADAMTS1; KIAA1346; METH1; A disintegrin and metalloproteinase with thrombospondin motifs 1; ADAM-TS 1; ADAM-TS1; ADAMTS-1; METH-1
Human Swiss Prot No.:	Q9UHI8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ADAMTS-1 Polyclonal Antibody detects endogenous levels of ADAMTS-1 protein.
Storage and Stability:	-20°C/1 year