

Antibodies Product Datasheet

Product Name:	ADAMTS-12 Antibody
Catalog Number:	C14304
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ADAMTS-12 . at AA range: 1100-1180
Concentration:	1 mg/ml
Gene Name:	ADAMTS12
Protein Name:	A disintegrin and metalloproteinase with thrombospondin motifs 12
Other Name:	ADAMTS12; A disintegrin and metalloproteinase with thrombospondin motifs 12; ADAM-TS 12; ADAM-TS12; ADAMTS-12
Human Swiss Prot No.:	P58397
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ADAMTS-12 Polyclonal Antibody detects endogenous levels of ADAMTS-12 protein.
Storage and Stability:	-20°C/1 year