

Antibodies Product Datasheet

Product Name:	ADAMTS-16 Antibody
Catalog Number:	C14308
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ADAMTS-16 . at AA range: 440-520
Concentration:	1 mg/ml
Gene Name:	ADAMTS16
Protein Name:	A disintegrin and metalloproteinase with thrombospondin motifs 16
Other Name:	ADAMTS16; KIAA2029; A disintegrin and metalloproteinase with thrombospondin motifs 16; ADAM-TS 16; ADAM-TS16; ADAMTS-16
Human Swiss Prot No.:	Q8TE57
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ADAMTS-16 Polyclonal Antibody detects endogenous levels of ADAMTS-16 protein.
Storage and Stability:	-20°C/1 year