

Antibodies Product Datasheet

Product Name:	ADAMTS-18 Antibody
Catalog Number:	C14310
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ADAMTS-18 . at AA range: 1030-1110
Concentration:	1 mg/ml
Gene Name:	ADAMTS18
Protein Name:	A disintegrin and metalloproteinase with thrombospondin motifs 18
Other Name:	ADAMTS18; ADAMTS21; A disintegrin and metalloproteinase with thrombospondin motifs 18; ADAM-TS 18; ADAM-TS18; ADAMTS-18
Human Swiss Prot No.:	Q8TE60
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ADAMTS-18 Polyclonal Antibody detects endogenous levels of ADAMTS-18 protein.
Storage and Stability:	-20°C/1 year