

Antibodies Product Datasheet

Product Name:	ADAMTS-19 Antibody
Catalog Number:	C14311
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ADAMTS19. AA range:251-300
Concentration:	1 mg/ml
Gene Name:	ADAMTS19
Protein Name:	A disintegrin and metalloproteinase with thrombospondin motifs 19
Other Name:	ADAMTS19; A disintegrin and metalloproteinase with thrombospondin motifs 19; ADAM-TS 19; ADAM-TS19; ADAMTS-19
Human Swiss Prot No.:	Q8TE59
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ADAMTS-19 Polyclonal Antibody detects endogenous levels of ADAMTS-19 protein.
Storage and Stability:	-20°C/1 year