

Antibodies Product Datasheet

Product Name:	ADAMTS-7 Antibody
Catalog Number:	C14318
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ADAMTS-7 . at AA range: 150-230
Concentration:	1 mg/ml
Gene Name:	ADAMTS7
Protein Name:	A disintegrin and metalloproteinase with thrombospondin motifs 7
Other Name:	ADAMTS7; A disintegrin and metalloproteinase with thrombospondin motifs 7; ADAM-TS 7; ADAM-TS7; ADAMTS-7; COMPase
Human Swiss Prot No.:	Q9UKP4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ADAMTS-7 Polyclonal Antibody detects endogenous levels of ADAMTS-7 protein.
Storage and Stability:	-20°C/1 year