

Antibodies Product Datasheet

Product Name:	ANT3 Antibody
Catalog Number:	C14334
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SLC25A6. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	SLC25A6
Protein Name:	ADP/ATP translocase 3
Other Name:	SLC25A6; ANT3; CDABP0051; ADP/ATP translocase 3; ADP; ATP carrier protein 3; ADP,ATP carrier protein, isoform T2; ANT 2; Adenine nucleotide translocator 3; ANT 3; Solute carrier family 25 member 6
Human Swiss Prot No.:	P12236
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ANT3 Polyclonal Antibody detects endogenous levels of ANT3 protein.
Storage and Stability:	-20°C/1 year