

Antibodies Product Datasheet

Product Name:	ANT4 Antibody
Catalog Number:	C14335
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SLC25A31. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	SLC25A31
Protein Name:	ADP/ATP translocase 4
Other Name:	SLC25A31; AAC4; ANT4; SFEC; ADP/ATP translocase 4; ADP; ATP carrier protein 4; Adenine nucleotide translocator 4; ANT 4; Solute carrier family 25 member 31; Sperm flagellar energy carrier protein
Human Swiss Prot No.:	Q9H0C2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ANT4 Polyclonal Antibody detects endogenous levels of ANT4 protein.
Storage and Stability:	-20°C/1 year