

Antibodies Product Datasheet

Product Name:	BST-1 Antibody
Catalog Number:	C14339
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BST1. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	BST1
Protein Name:	ADP-ribosyl cyclase 2
Other Name:	BST1; ADP-ribosyl cyclase 2; Bone marrow stromal antigen 1; BST-1; Cyclic ADP-ribose hydrolase 2; cADPr hydrolase 2; CD antigen CD157
Human Swiss Prot No.:	Q10588
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BST-1 Polyclonal Antibody detects endogenous levels of BST-1 protein.
Storage and Stability:	-20°C/1 year