

Antibodies Product Datasheet

Product Name:	ARF GAP3 Antibody
Catalog Number:	C14346
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ARF GAP3 . at AA range: 280-360
Concentration:	1 mg/ml
Gene Name:	ARFGAP3
Protein Name:	ADP-ribosylation factor GTPase-activating protein 3
Other Name:	ARFGAP3; ARFGAP1; ADP-ribosylation factor GTPase-activating protein 3; ARF GAP 3
Human Swiss Prot No.:	Q9NP61
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ARF GAP3 Polyclonal Antibody detects endogenous levels of ARF GAP3 protein.
Storage and Stability:	-20°C/1 year