

Antibodies Product Datasheet

Product Name:	GPT2 Antibody
Catalog Number:	C14370
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from GPT2 . at AA range: 190-270
Concentration:	1 mg/ml
Gene Name:	GPT2
Protein Name:	Alanine aminotransferase 2
Other Name:	GPT2; AAT2; ALT2; Alanine aminotransferase 2; ALT2; Glutamate pyruvate transaminase 2; GPT 2; Glutamic--alanine transaminase 2; Glutamic--pyruvic transaminase 2
Human Swiss Prot No.:	Q8TD30
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPT2 Polyclonal Antibody detects endogenous levels of GPT2 protein.
Storage and Stability:	-20°C/1 year