

Antibodies Product Datasheet

Product Name:	ADH7 Antibody
Catalog Number:	C14383
Reactivity:	Human, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ADH7. AA range:211-260
Concentration:	1 mg/ml
Gene Name:	ADH7
Protein Name:	Alcohol dehydrogenase class 4 mu/sigma chain
Other Name:	ADH7; Alcohol dehydrogenase class 4 mu/sigma chain; Alcohol dehydrogenase class IV mu/sigma chain; Gastric alcohol dehydrogenase; Retinol dehydrogenase
Human Swiss Prot No.:	P40394
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ADH7 Polyclonal Antibody detects endogenous levels of ADH7 protein.
Storage and Stability:	-20°C/1 year