

Antibodies Product Datasheet

Product Name:	DD1 Antibody
Catalog Number:	C14393
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from DD1 . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	AKR1C1
Protein Name:	Aldo-keto reductase family 1 member C1
Other Name:	AKR1C1; DDH; DDH1; Aldo-keto reductase family 1 member C1; 20-alpha-hydroxysteroid dehydrogenase; 20-alpha-HSD; Chlordecone reductase homolog HAKRC; Dihydrodiol dehydrogenase 1/2; DD1/DD2; High-affinity hepatic bile acid-binding protein; HB
Human Swiss Prot No.:	Q04828
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DD1 Polyclonal Antibody detects endogenous levels of DD1 protein.
Storage and Stability:	-20°C/1 year