

Antibodies Product Datasheet

Product Name:	DD2 Antibody
Catalog Number:	C14394
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AKR1C2. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	AKR1C2
Protein Name:	Aldo-keto reductase family 1 member C2
Other Name:	AKR1C2; DDH2; Aldo-keto reductase family 1 member C2; 3-alpha-HSD3; Chlordecone reductase homolog HAKRD; Dihydrodiol dehydrogenase 2; DD-2; DD2; Dihydrodiol dehydrogenase/bile acid-binding protein; DD/BABP; Trans-1; 2-dihydrobenzene-1,2-diol
Human Swiss Prot No.:	P52895
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DD2 Polyclonal Antibody detects endogenous levels of DD2 protein.
Storage and Stability:	-20°C/1 year