

Antibodies Product Datasheet

Product Name:	ASAH3 Antibody
Catalog Number:	C14401
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ASAH3 . at AA range: 100-180
Concentration:	1 mg/ml
Gene Name:	ACER1
Protein Name:	Alkaline ceramidase 1
Other Name:	ACER1; ASAH3; Alkaline ceramidase 1; AlkCDase 1; Alkaline CDase 1; Acylsphingosine deacylase 3; N-acylsphingosine amidohydrolase 3
Human Swiss Prot No.:	Q8TDN7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ASAH3 Polyclonal Antibody detects endogenous levels of ASAH3 protein.
Storage and Stability:	-20°C/1 year