

Antibodies Product Datasheet

Product Name:	ASAH3L Antibody
Catalog Number:	C14402
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ASAH3L . at AA range: 50-130
Concentration:	1 mg/ml
Gene Name:	ACER2
Protein Name:	Alkaline ceramidase 2
Other Name:	ACER2; ASAH3L; PP11646; Alkaline ceramidase 2; AlkCDase 2; Alkaline CDase 2; haCER2; Acylsphingosine deacylase 3-like; N-acylsphingosine amidohydrolase 3-like
Human Swiss Prot No.:	Q5QJU3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ASAH3L Polyclonal Antibody detects endogenous levels of ASAH3L protein.
Storage and Stability:	-20°C/1 year