

Antibodies Product Datasheet

Product Name:	APHC Antibody
Catalog Number:	C14411
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PHCA. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	ACER3
Protein Name:	Alkaline ceramidase 3
Other Name:	ACER3; APHC; PHCA; Alkaline ceramidase 3; AlkCDase 3; Alkaline CDase 3; Alkaline dihydroceramidase SB89; Alkaline phytoceramidase; aPHC
Human Swiss Prot No.:	Q9NUN7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	APHC Polyclonal Antibody detects endogenous levels of APHC protein.
Storage and Stability:	-20°C/1 year