

Antibodies Product Datasheet

Product Name:	A1BG Antibody
Catalog Number:	C14418
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from A1BG . at AA range: 100-180
Concentration:	1 mg/ml
Gene Name:	A1BG
Protein Name:	Alpha-1B-glycoprotein
Other Name:	A1BG; Alpha-1B-glycoprotein; Alpha-1-B glycoprotein
Human Swiss Prot No.:	P04217
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	A1BG Polyclonal Antibody detects endogenous levels of A1BG protein.
Storage and Stability:	-20°C/1 year