

Antibodies Product Datasheet

Product Name:	ENaC δ Antibody
Catalog Number:	C14438
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SCNN1D. AA range:411-460
Concentration:	1 mg/ml
Gene Name:	SCNN1D
Protein Name:	Amiloride-sensitive sodium channel subunit delta
Other Name:	SCNN1D; DNACH; Amiloride-sensitive sodium channel subunit delta; Delta-NaCH; Epithelial Na(+) channel subunit delta; Delta-ENaC; ENaCD; Nonvoltage-gated sodium channel 1 subunit delta; SCNED
Human Swiss Prot No.:	P51172
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ENaC δ Polyclonal Antibody detects endogenous levels of ENaC δ protein.
Storage and Stability:	-20°C/1 year