

Antibodies Product Datasheet

Product Name:	AMPD1 Antibody
Catalog Number:	C14443
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AMPD1. AA range:261-310
Concentration:	1 mg/ml
Gene Name:	AMPD1
Protein Name:	AMP deaminase 1
Other Name:	AMPD1; AMP deaminase 1; AMP deaminase isoform M; Myoadenylate deaminase
Human Swiss Prot No.:	P23109
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AMPD1 Polyclonal Antibody detects endogenous levels of AMPD1 protein.
Storage and Stability:	-20°C/1 year